INSTRUCTIUNI DE UTILIZARE INSTRUCTIONS FOR USE

Brand: SIFRA Producator / Manufacturer: MERCAN OPTIK ITHALAT IHRACAT SAN.
VE TIC. LTD. STI.

Clasa / Class: I (MDR 2017/745) GMDN: 32816 Ref: ISO 12870

■■ ROMANA

■■ ENGLISH

1. DOMENIU DE UTILIZARE

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1. INTENDED USE

SIFRA spectacle frames are Class I medical devices (Rule 1, Chapter III, Annex VIII of Regulation (EU) 745/2017).

They are intended exclusively for use with prescription ophthalmic lenses, to support ocular prostheses and correct visual defects.

The product complies with Regulation (EU) 2017/745 on medical devices and with standard EN ISO 12870 (ophthalmic frames — requirements and test methods).

Not intended for use as sunglasses, protective eyewear, or any other non-medical purpose.

2. WARNINGS

The frame is delivered with demonstration lenses, which are not intended for use and serve only to maintain the shape and allow fitting under realistic conditions.

This is not personal protective equipment — do not use it as such.

Keep out of reach of children.

If visible damage is present (cracks, deformation), discontinue use and consult an authorised optician.

3. CLEANING AND MAINTENANCE

Store the frames at a temperature not exceeding 90°C, in the provided case.

Avoid prolonged exposure to sunlight in enclosed spaces (e.g. car dashboard).

If in contact with salt water, pool water, cosmetics or hairspray, wipe immediately with a soft damp cloth.

Use only products recommended by your optician: microfibre cloth, pH-neutral solutions or special lens cleaning wipes.

Avoid solvents (acetone, alcohol) and chemical detergents that may damage the frame materials.

In case of wear (scratches, stains), replace components only with original parts, at authorised optical centres.

4. CE MARKING

The CE marking affixed to the product certifies compliance with the essential requirements of Regulation (EU) 2017/745 on medical devices.

The model identification code and product dimensions are found on the inner side of the temples.

Basic UDI-DI: 8680277BXPLSGC01XZ / 8680277BXMTLGC01X9 / 8680277BXMIXGC01WW

5. MANUFACTURER

MERCAN OPTIK ITHALAT IHRACAT SAN. VE TIC. LTD. STI.

Ikiteilli O.S.B. Mutsan Sitesi M4 Blok No:10

Basaksehir — Istanbul / Turkey

SRN: TR-MF-000052576

6. COMPLAINTS AND REPORTS

Romanian importer: Sifra Vision SRL

Str. G-ral Alexandru Cernat nr. 188, Galati, Romania

E-mail: office@sifravision.ro | Web: www.sifravision.ro

Any serious incident related to the use of this medical device must be reported to the national competent authority (ANMDMR).